



## 2.1 mm ID Micro Carbamates Column

Ideal for  
LC-MS  
Confirmation

Post-column derivatization of carbamates with hydroxide and OPA is the method of choice for the detection of carbamates because of its high selectivity and high sensitivity (US EPA 531.1 & AOAC 985.23).

To confirm the finding of carbamates in the tested samples, there are various possibilities:

- 1) Alternate columns to give different retention time or order of elution.
- 2) Fluorescence scan of the isoindole end-products.
- 3) Mass-spectrum of the carbamates.

Options 1 & 2 are already available to many Pickering customers. For option 3, the resolution of the carbamates has always been an issue with the necessary smaller ID column. The new 2.1 mm ID Pickering Carbamates Column provides the resolution. Because the new column contains the same durable phase as our regular Carbamates Column, it will have minimum leaching (critical for the mass spectrometer).


### METHOD

#### Analytical Conditions

COLUMN: 2.1 x 100 mm Carbamates Column  
Catalog No. 1821100  
TEMPERATURE: 42 °C  
FLOW RATE: 0.2 mL/min  
MOBILE PHASE: A=H<sub>2</sub>O; B=MeOH

#### Post-column Conditions

POST-COLUMN SYSTEM: Pinnacle PCX or Vector PCX  
REAGENT 1: CB 130  
REACTOR 1: 100 µL, 100 °C  
REAGENT 2: OPA, Thiofluor in CB 910  
REACTOR 2: 20 µL, Ambient  
FLOW RATE: 60 µL/min  
DETECTION: Fluorometer  
λ<sub>ex</sub>: 330 nm  
λ<sub>em</sub>: 464 nm



Time [min]	A [%]	B [%]
0.0	100	0
0.1	88	12
15.0	60	40
20.0	55.5	44.5
27	20	80
27.1	0	100
33	0	100
45.1	88	12