

Alprostadil (Prostaglandin E1) is a drug that has vasodilation properties and is used to treat erectile dysfunction (ED) and other medical conditions. This medication is available in injectable form and in suppository form.

Since Prostaglandin E1 has low UV absorbance, the analysis of the formulations could be difficult. The addition of post-column derivatization increases the sensitivity and specificity of the analysis.

## Method

### Analytical Conditions

**Analytical Column:** Reversed-phase  $C_{18}$  column (150 x 4.6 mm)

**Temperature:** 37 °C

**Flow Rate:** 1 mL/min

**Mobile Phase:** 30% Acetonitrile – 70% of 0.02 M Potassium Phosphate monobasic (adjusted to pH 3)

Or alternative conditions:

**Mobile Phase:** 25% Acetonitrile – 75% of 0.0067 M Potassium Phosphate buffer pH 6.3

**Injection Volume:** 20  $\mu$ L

### Post-Column Conditions

**Post-Column System:** Onyx PCX or Pinnacle PCX

**Reactor Volume:** 2 mL

**Reactor Temperature:** 60 °C

**Reagent:** 1 mol/L Potassium Hydroxide

**Reagent Flow Rate:** 1 mL/min

**Detection:** UV 278 nm

### Calibration

**Alprostadil:** 0.1  $\mu$ g/mL – 10  $\mu$ g/mL

**$\beta$ -Naphthol (Internal Standard):** 0.25  $\mu$ g/mL – 10  $\mu$ g/mL

Make stock solutions of Alprostadil and  $\beta$ -Naphthol in anhydrous Ethanol. Prepare working standards by appropriate dilution of stock solution with mobile phase.




Fig 1. Calibration curve for Alprostadil




Fig 2. Calibration curve for  $\beta$ -Naphthol




Fig 3. Chromatogram of Alprostadil and  $\beta$ -Naphthol (Internal Standard)