

Voglibose is an Alpha-Glucosidase inhibitor widely used for the treatment of diabetes. Alpha-glucosidase inhibitors are agents that delay the glucose absorption at the intestinal level and thereby prevent a sudden surge of glucose after a meal. Voglibose is the safest and most effective drug of its class.

Since Voglibose has no UV chromophore, post-column derivatization is employed to produce a fluorescent derivative.

This abstract describes a very sensitive and robust analytical method for the analysis of Voglibose in pharmaceutical tablets. Simple sample preparation and fast analysis time allow for using this method in high throughput environments.

### Method

#### Analytical Conditions

**Column:** Restek Pinnacle II Amino,  
5  $\mu$ m, (250 x 4.6) mm,  
Catalog # 9217575

**Temperature:** 35 °C

**Flow Rate:** 0.6 mL/min

**Mobile Phase:** Sodium phosphate buffer, 20 mM  
pH 6.5 / Acetonitrile (37 : 63)

**Injection Volume:** 50  $\mu$ L

**Note:** We strongly recommend that the column be flushed with acetonitrile water (80:20) for twenty minutes before making any injections.

#### Sample Preparation

Crush 5 tablets and mix with 25 mL of mobile phase. Sonicate for 10 min and filter liquid portion through 0.45  $\mu$ m filter. Put in HPLC autosampler vial and inject 50  $\mu$ L.

#### Post-Column Conditions

**Post-Column System:** Onyx PCX, Pinnacle PCX or Vector PCX

**Heated Reactor Volume:** 3.5 mL

**Temperature:** 100 °C

**Cooling Coil:** 0.15 mL (at room temperature)


**Reagent:** Taurine (6.25 g), Sodium Periodate (2.56 g)  
in 1000 mL of water

**Flow Rate:** 0.6 mL/min


**Detection:** FLD  
 $\lambda_{\text{ex}}$ : 350 nm,  $\lambda_{\text{em}}$ : 430 nm

Repeatability Studies for  
Different Concentration Levels


	0.5 ppm	100 ppm
Average RT, min	21.25	21.26
RSD, %, N = 6	0.36	0.08
Average Peak Area	9.22	1,562.69
RSD, %, N = 6	1.48	0.79




Calibration Curve for Analytical Range 0.5-200 ppm



Chromatogram of Voglibose Standard, 50 ppm, 50 µL Injection



Chromatogram of Voglibose Standard, 0.5 ppm, 50 µL Injection



Chromatogram of Voglibose Tablets (Volix™, 0.2 mg), 50 µL Injection

