



**Restrictor, 6 cm, 200 psi
PN 1100-0205**

Dear valued customer,

After recent research demonstrating enhanced performance capabilities, we are now offering a replacement for part # 1100-0205 (OPA Flow Restrictor, 300 psi). Effective immediately, we will be discontinuing 1100-0205 and replacing with part # 1100-0215 (Restrictor, PEEK Tubing, 200 psi). The replacement PEEK Tubing Flow Restrictor is an improvement to our old flow restrictor and comes with the following benefits:

- More robust construction material
- Longer usable life that translates to more sample runs
- Fewer connections creating less points of failure
- Less expensive price point

Please contact Pickering Laboratories Support at (650) 694-6700 or support@pickeringlabs.com if you have any questions.

Obsolete Part #	Replacement Part #	Replacement Description
1100-0205	1100-0215	Restrictor, PEEK Tubing, 200 psi

Regards,

David Mazawa
Technical Support Chemist
Pickering Laboratories, Inc.
1280 Space Park Way
Mountain View, CA 94043 USA
Phone: (650) 694-6706