



## TEST MIXTURE Specifications & Handling Information

**Catalog No.** 1700-0063 (5.0 mL)

**Description:** Carbamate Test Mixture, in Methanol, 2.5 µg/mL

---

### ▪ USAGE

Test Mixtures are the qualitative standards and intended to be used to establish an elution profile and for troubleshooting. The Carbamate Test Mixture is a methanol solution of eleven Carbamate insecticides and 1-Naphthol at a concentration of approximately 2.5 µg/mL.

### ▪ STORAGE

Store frozen.

### ▪ CONSTITUENTS

Aldicarb  
Aldicarb Sulfone  
Aldicarb Sulfoxide  
BDMC (4-Bromo-3,5-dimethylphenyl-N-methylcarbamate)  
Carbaryl  
Carbofuran  
3-Hydroxycarbofuran  
Methiocarb  
Methomyl  
1-Naphthol  
Oxamyl  
Propoxur