



TEST MIXTURE
Specifications & Handling Information

Catalog No. **1700-0080 (5.0 mL)**

Description: Glyphosate Test Mixture, in Water, 2.5 µg/mL

▪ **USAGE**

Test Mixtures are the qualitative standards and intended to be used to establish an elution profile and for troubleshooting. The Glyphosate Test Mixture is a water solution of Glyphosate and Aminomethylphosphonic Acid (AMPA) at a concentration of approximately 2.5 µg/mL.

▪ **STORAGE**

Store frozen.

▪ **CONSTITUENTS**

Glyphosate
AMPA (Aminomethyl-phosphonic acid)