



Specifications & Storage Information

Catalog No. **1700-0190** (5 x 1 mL)

Description: Carbamates Calibration Standard in Acetonitrile

CONSTITUENTS

Aldicarb Sulfone 100 ug/mL
Aldicarb Sulfoxide 100 ug/mL
Oxamyl 100 ug/mL
Methomyl 100 ug/mL
3-Hydroxycarbofuran 100 ug/mL
Aldicarb 100 ug/mL
Propoxur 100 ug/mL
Carbofuran 100 ug/mL
1-Naphthol 100 ug/mL
Carbaryl 100 ug/mL
Methiocarb 100 ug/mL

Catalog No. **1700-0185** (5 x 1 mL)

Description: Carbamates Internal Standard (BDMC) in Acetonitrile

CONSTITUENTS

4-Bromo-3,5-dimethylphenyl-N-methylcarbamate (BDMC) 100 ug/mL

STORAGE

Frozen at -15°C.

[DOWNLOAD CERTIFICATES OF ANALYSIS](https://www.pickeringlabs.com/library/certificate-of-analysis/)

<https://www.pickeringlabs.com/library/certificate-of-analysis/>